



Finger Injury Kit

Model : 740-640

Every finger injury is different and requires specific treatment. Please consult your physician prior to starting therapy.

Designed to help relieve pain and reduce swelling caused by sprains, strains, arthritis, and more.

Cold Therapy:

For maximum performance, place Therm-O-Gel in freezer for at least 2 hours prior to use. Always return to freezer following treatment.

Suggestions:

Use Therm-O-Gel 3-4 times a day for 10 minutes per treatment until swelling and soreness subsides. Slide gauze over finger before applying the gel pack. Place the unit around your finger and use the strap to fasten it, ensuring that there isn't any direct pressure on the wound. Strap can be removed and repositioned if necessary. If redness, pain, or swelling persist, consult a physician immediately.

www.proactive-painmanagement.com

Manufactured and printed in China for /  A.M.G. Medical Inc. 8505 Dallon,
Fabriqué et imprimé en Chine pour : Montréal, QC H4T 1V5 Canada • 1-800-363-2381

Distributed in the USA by / **AMG Medical, Inc.** (USA) 8396 State Route 9
Distribué aux É.-U. par : West Chazy, NY 12992 USA • 1-888-412-4992



Trousse pour blessure au doigt

Modèle : 740-640

Chaque blessure aux doigts est différente et nécessite un traitement spécifique. Veuillez consulter votre médecin avant de commencer le traitement.

Conçu pour aider à soulager la douleur et à réduire l'enflure causée par les entorses, les foulures, l'arthrite, etc.

Thérapie par le froid:

Pour des performances maximales, placez Therm-O-Gel au congélateur pendant au moins 2 heures avant de l'utiliser. Remettez toujours la compresse au congélateur après le traitement.

Suggestions:

Utilisez Therm-O-Gel 3 à 4 fois par jour pendant 10 minutes jusqu'à ce que l'enflure et la douleur disparaissent. Enfilez la gaze sur le doigt et insérez celui-ci dans la compresse. Utilisez la sangle pour la fixer, en veillant à ce qu'il n'y ait pas de pression directe sur la blessure. La sangle peut être retirée et repositionnée au besoin. Si les rougeurs, douleurs ou gonflements persistent, consultez immédiatement un médecin.

www.proactive-painmanagement.com