

# 770-300 MEDPRO AQUACARE COMMODE WITH INFECTION CONTROL FRIENDLY SEAT & ARM PADS



**770-300**  
**MEDPRO AQUACARE COMMODE WITH**  
**INFECTION CONTROL FRIENDLY SEAT & ARM PADS**

Item #	AMG Part #	Item Description	Qty Per 770-300	Pkg Qty
1.	<b>B01-259</b>	Blue Backrest Pad with Hardware (4 Brackets & 8 Philips Screws)	1 EA	1 EA
2.	<b>B01-451</b>	Off-White Arm Pad with Hardware (Off-White HDPE Pad & 2 Philips Screws M5x32mm) - 2ea/set	1 ST	1 ST
3.	<b>B01-450</b>	Complete Armrest Set for MedPro Aquacare Commode (2 Off-White Arm Pads & 2 Philips Screws M5x32mm, Chrome Plated Arm Bar Left & Right, 2 male grey holders & 4 Setscrews)	1 PR	1 PR
4.	<b>B01-481</b>	Male Hardware set for Armrest (2 male grey holders & 4 Setscrews)	1 ST	1 ST
5.	<b>B01-482</b>	Female Hardware set for Armrest (2 female grey holders, 2 Philips Screws M6x30mm & 2 Locknuts)	1 ST	1 ST
6.	<b>B01-171</b>	Grey Plastic Caster with Locking System (4ea / set)	1 ST	1 ST
7.	<b>A01-411</b>	Grey Tightening Knob (2ea / set)	1 ST	1 ST
8.	<b>A01-436</b>	Blue Soft Seat Cover	1 EA	1 EA
9.	<b>B01-261</b>	Off-White Blow-Molded Seat with Pail Support Brackets	1 EA	1 EA
*9.	<b>770-380</b>	Off-White Blow-Molded Seat with Bed Pan Support Brackets <b>(Accessory)</b>	1 EA	1 EA
10.	<b>B01-710</b>	Pail and Cover	1 EA	1 EA
11.	<b>A01-441</b>	Grey Plastic Foot-Support	1 EA	1 EA