


770-375 EURO COMMODE (WITH 4 LOCKING CASTERS)



Item #	AMG Part #	Item Description	Qty Per Unit	Pkg Qty
1.	B01-231	Armrest with Hardware	1 PR	1 PR
2.	B01-710	Pail and Cover	1 EA	1 EA
3.	A01-412	Tightening Knob for Backrest	2 EA	1 EA
4.	B01-260	Toilet Seat with Pail Bracket	1 EA	1 EA
5.	770-364	5" Locking Caster	4 EA	1 EA
6.	770-359	Bedpan Holder (Optional Accessory)	1 EA	1 EA
7.	A01-431	Padded Seat, Blue	1 EA	1 EA
* 8.	B02-306	Swing-Away Footrests, Tool-free Adjustable (Left & Right)	1 PR	1 PR
9.	A01-432	Velcro Heel Strap	1 EA	1 EA

** For more details see B02-306 file*