


KD BASE I.V. STAND
775-840
STAINLESS STEEL POLE
(WITHOUT ACCESSORIES)



KD BASE I.V. STAND
775-840
STAINLESS STEEL POLE
(WITHOUT ACCESSORIES)

Item #	AMG Part #	Item Description	Qty Per Unit	Sold by
1.	A03-560	4 Hooks Top (Chromed Steel)	1 EA	1 EA
2.	A03-460	Orange Tightening Knob - 3ea / set	1 ST	1 ST
3.	B03-300	Orange Telescoping Mechanism with Orange Knob	1 ST	1 ST
4.	B03-002	Stainless Steel KD Pole (lower+middle+upper pole)	1 BX	1 BX
5.	B03-220	Patient Handle (Chromed Steel) with Orange Knob – Sold like Accessory	1 BX	1 BX
6.	B03-230	Drainage Bag Hook (Chromed Steel) with Orange Knob – Sold like Accessory	1 BX	1 BX
7.	B03-071	Chromed-Finish KD Base Legs with Hardware (5 Chromed-Finish Legs, 5 pop-in 3" Grey Casters, 1 Stainless Steel Cover, 5 Screws M8x20mm, 5 Lock Washers for M8, 1 Bolt M10x50mm, 1 Washer for M10 & 1 Rubber Ring)	1 BX	1 BX
8.	B03-110	Pop-in 3" Grey Caster	1 EA	1 EA
9.	B03-870	Hardware for KD Base (1 Stainless-Steel Cover, 5 Screws M8x20mm, 5 Lock Washers for M8, 1 Bolt M10x50mm, 1 Washer for M10)	1 BX	1 BX